



# Clinical Protocol Vital Pulp



# Clinical Protocol Necrotic, Apical Periodontitis



# Clinical Protocol Retreatment

Minimal apical box  
# 35/40/0.02/0.04 taper  
(curvature)

1 % NaOCl

Ultrasonication 30-60s

17 % EDTA or 10 %- 20 %  
Citric Acid

Ultrasonication 30-60s

Drying with paper points

Clorhexidin digluconat 2 %

Ultrasonication 30-60s

Obturation

IC Medication 3-4 weeks



**TECHNOMEDICS NORGE A.S**  
DENTAL & KIRURGI