



# Get to the Point.

## Elements Diagnostic Unit and Apex Locator

When you're doing a root canal, you need to get to the point - the apex - predictably, efficiently and accurately. SybronEndo's **Elements Diagnostic Unit and Apex Locator** provides the convenience and accuracy of a high quality pulp tester and apex locator in one system. What's more, the satellite display conveniently brings critical data into your field of vision for a faster more stable reading.





	PART #	DESCRIPTION
A	973-0310	File Clip
B	973-0307	Bifurcated Probe
C	973-0309	Crown Test Tip (Diagnostic Unit only)
D	973-0308	Validity Probe (Diagnostic Unit Only)
E	973-0313	Lip Hook (5-pack)
F	973-0311	Elements Patient Lead Cord
G	973-0303	Satellite with Satellite Cord
H	973-0300	Elements Diagnostic Unit
H	973-0301	Elements Apex Locator

FEATURES	BENEFITS
<ul style="list-style-type: none"> <li>Versatile Satellite display</li> </ul>	<ul style="list-style-type: none"> <li>Quick, easy reading and unparalleled mounting flexibility - on patient bib, instrument tray, or microscope</li> </ul>
<ul style="list-style-type: none"> <li>Auto-power off</li> </ul>	<ul style="list-style-type: none"> <li>Conserves battery life</li> </ul>
<ul style="list-style-type: none"> <li>Full range audio speakers</li> </ul>	<ul style="list-style-type: none"> <li>Eliminates annoying beeps</li> </ul>
<ul style="list-style-type: none"> <li>Non-glare tilt screen with digital and graphic feedback</li> </ul>	<ul style="list-style-type: none"> <li>Large and easy to read</li> </ul>
<ul style="list-style-type: none"> <li>Autoclavable medical-grade cords and connectors</li> </ul>	<ul style="list-style-type: none"> <li>Durable and easy to clean</li> </ul>
<ul style="list-style-type: none"> <li>Wide array of gold plated probes</li> </ul>	<ul style="list-style-type: none"> <li>Resists micro-corrosion for more accurate feedback</li> </ul>
<ul style="list-style-type: none"> <li>Automatic calibration</li> </ul>	<ul style="list-style-type: none"> <li>Accurate analysis in the presence of EDTA, sodium hypochlorite, blood, saline, water and hydrogen peroxide</li> </ul>
<ul style="list-style-type: none"> <li>Validity Scanner</li> </ul>	<ul style="list-style-type: none"> <li>Pain free electric pulp testing of teeth regardless of shape, size or abnormalities</li> </ul>

**SybronEndo** 1717 West Collins Avenue, Orange, CA 92867 Tel 800.346.ENDO Tel 714.516.7979 Fax 714.516.7936 sybronendo.com

**SybronEndo Europe** Basicweg 20, 3821 BR Amersfoort, The Netherlands Tel +31 33 453 61 80 Fax +31 33 453 61 53 customer.endo@sybrondental.com

**SybronEndo Australia** Unit 10/112-118 Talavera Road, North Ryde NSW 2113 Tel 1800 052 030 Fax 02 8875 8175 sybronendoaustralia@sybrondental.com